

BIOBASE

Potassium (K⁺) Test Kit (Enzymatic)

【NAME】

Name: Potassium (K) Test Kit (Enzymatic)

【Package】

R1:60mL×2 R2:20mL×2 R1:18mL×4 R2:6mL×4
R1:18mL×1 R2:6mL×1 R1:60mL×1 R2:20mL×1
R1:40mL×3 R2:20mL×2 R1:48mL×2 R2:16mL×2
R1:45mL×2 R2:15mL×2 R1:90mL×1 R2:30mL×1
2×200Tests

【INTENDE USE】

The reagent is intended for the in vitro quantitative determination of Potassium (K) in human serum, plasma on both manual and automatic systems.

(1) Serum Potassium increase often occur in the following situations:


① Acute renal failure, potassium excretion dysfunction when urine intoxication; ② the adrenal cortex hypofunction; ③ a variety of reasons of respiratory acidosis and metabolic acidosis; ④ Severe hemolysis, a large number of input banked blood, crush syndrome, burn, etc.; ⑤ Extensive use of potassium drugs; ⑥ High permeability dehydration.

(2) Serum Potassium reduce often appears in the following situations:

① such as severe infection, chronic wasting disease long-term inappetence and fast after the surgery time is too long and did not pay attention to the potassium supplement; ② Adrenal cortex hyperfunction. long-term heavy use of adrenal cortex hormone; Many long-term use of diuretics; acute renal failure stage by anuresis into the polyuria; ③ Alkali poisoning; Diabetic patients using insulin or insulin and glucose as energy mixture is used; ④ A large number of input without liquid potassium, plasma dilution, lower the serum potassium.

【METHODOLOGY】

Potassium detection by potassium dependent pyruvate kinase catalytic enol phosphate type pyruvate (PEP) enzyme kinetics reaction. Its product (pyruvic acid salt) and NADH reaction generated NAD⁺ under the action of lactate dehydrogenase (LDH), the product (NAD⁺) absorbance value is proportional to the concentration of potassium in the 340 nm.



【REAGENT COMPOSITION】

Reagent 1 (R1) :

Tris buffer PH8.2	250mmol/L
Cryptand	12mmol/L
PEP	≥3.0mmol/L
ADP	≥3.2mmol/L
α-ketone glutaric acid	≥1.2mmol/L
NADH	≥0.35mmol/L
GLDH	≥11U/mL
PK	≥1.2U/mL

Reagent 2 (R2) :

LDH ≥65U/mL

【STABILITY AND STORAGE】

Unopened, avoid light preservation in 2~8℃, valid for 12 months;

Opened, avoid contamination preservation in 2~8℃, valid for 1 month.

【SPECIMEN COLLECTION AND HANDLING】

Serum or lithium heparin anticoagulant blood plasma.

Serum stability: 2~8℃ preservation stability in 3 days;

When the bilirubin concentration of sample ≤665μmol/L, Hb ≤1g/L, TG ≤24.2mmol/L, was not observed clearly disturbance.

【APPLICABLE INSTRUMENT】

BIOBASE、HITACHI、OLYMPUS、BECKMAN、GLAMOUR、ABBOTT、TOSHIBA, etc. fully automatic biochemical analyzer.

【TESTING SPECIFICATION】

1. Use requirement: Liquid reagent use after opening

2. Test condition:

Wavelength	340nm
Cuvette light path	1.0cm
Temperature	37℃
Assay Type	fix-time

3. Calibration procedure:

Two calibration with deionized water and K calibration. Calibration another purchase, it is recommended to use RANDOX, or other traceability calibration.


4. Quality control program:

Determination of K quality control products, test results within the scope of quality control can test the samples.

5. OPERATION STEPS

Double reagent operation:

R1: reagent 1 R2: reagent 2 S: Calibrator U: Sample



Distilled water、U or S:50μl; R2:600μl Recording absorbance A1 Recording absorbance A2
value, calculate ΔA/min

【CALCULATION】

Sample ΔA/min
K sample concentration = $\frac{\text{Sample } \Delta A / \text{min}}{\text{Calibrator } \Delta A / \text{min}}$ × Standard concentration

【REFERENCE RANGE】

3.5~5.1mmol/L (13.7~19.9mg/dL)

Unit conversion μmol/L × 3.9g/mol = mg/dL

By clinical trials, choose no less than 100 newborn or adults blood specimens, tested by automatic biochemical analyzer, and then processing the testing value with statistical method, calculating out the reference range.

【THE LIMITATIONS OF TESTING RESULTS】

K testing is just one of the standard that clinician diagnose the patient. Clinical physicians should according to patients' bodies, history and other diagnostic program, to get comprehensive judgment.

【THE INTERPRETATION OF TEST RESULTS】

Human error, the processing of specimen, analysis instrument deviation, etc. all can affect the measurement result; When one sample deviates from

the expected value too far, need to be tested again.

【PERFORMANCE INDEX】

1. Reagent blank absorbance ≥ 1.7000 (340nm, 1cm optical path).
2. Precision: repeatability $CV \leq 10\%$; batch variations $R \leq 10\%$.
3. Accuracy: relative deviation $\leq 10\%$.
4. Linearity range: 0~10mmol/L(7.8~39.0mg/dL), $r \geq 0.990$.
5. Stability: All package reagent, non corrosive gas and avoid light, preservation in 2~8 °C, stable 12 months.

【ATTENTION】

1. When the simultaneous determination of sodium and potassium project, please first test Na.
2. Different batches reagents cannot mix, when replacing reagents batch number, please calibration again.
3. The opened reagent should be airtight store according to the specified method, never use expired.
4. Liquid waste disposal: Suggest follow local regulations.

【Production enterprise licence number】: 鲁食药监械生产许 20080002 号

【Registration certificate number】: 鲁食药监械(准)字 2010 第 2400280 号

【Product standard number】: YZB/鲁 0206-2010

Biobase Bioindustry(Shandong)Co.,Ltd

Add: BIOBASE Biological industrial park, Mingshui Economic development zone, Jinan, China, 250200

Tel: 0086-0531-91219801 Fax: 0086-0531-81219804